

Introduction:

Dr. Jonathan Nelson

OHSU Dept. of Nephrology and Hypertension

Employment opportunity

# Basic gene cloning: plasmids and transformation

- Plasmid biology
- How to add DNA to a plasmid
- Transformation of *E. coli*

The key role of *E. coli* and an RNA expression plasmid in mRNA vaccine production (Pfizer/BioNTech, specifically)

<https://www.nytimes.com/interactive/2021/health/pfizer-coronavirus-vaccine.html?action=click&module=Spotlight&pgtype=Homepage>

-

# Readings:

- 1) *22 MC4 Plasmids*. History of plasmid development, plasmid replicons, care of E. coli and plasmids, specialized plasmids, transformation, alpha complementation.
- 2) *25 MC4 Antibiotics*. Summary of several antibiotics used for plasmid maintenance, and how they work
- 3) *24 MC4 Cloning in Plasmids*. Practical aspects of cloning and transformation, plasmid dephosphorylation. Discussion of recombinational and ligation independent cloning.
- 4) Gibson assembly 2009. Report demonstrating easy DNA assembly using overlapping sequences.
- 5) IVA (in vitro assembly) 2016. Cloning without ligation.
- 6) Link: [www.addgene.org](http://www.addgene.org), a nonprofit plasmid repository

**Cloning:** isolation and propagation of a specific piece of DNA

**Cloning vector:** a movable DNA element that is modified to contain another piece of DNA

**General features:**

- Easy to isolate and purify
- Can be manipulated in the test tube
- Can be moved into a new cell
- Selectable: easy to propagate and test for presence or absence

**Types of cloning vectors:**

- Plasmid
- Virus
- Artificial chromosome
- transposon

# Bacterial plasmid: Non-chromosomal DNA molecule


- Small, circular, supercoiled
- Replicates independently of the bacterial chromosome
- Copy number in cell is maintained through an origin of replication (replicon)
- Selectable through a gene that is beneficial for the host bacterium, eg. antibiotic resistance

# A useful plasmid

1. Is relatively small
  - Gives higher copy #, stability, and transforming efficiency
  - Can accept larger pieces of DNA
  - Easier to handle (less susceptible to breakage)
2. Has a selectable marker (e.g. antibiotic resistance)
3. Has single sites for many restriction enzymes

# Plasmid: basic parts list

gene that encodes  
antibiotic resistance  
(for plasmid  
maintenance)



Cloning site

Replicon/origin  
of replication  
(for plasmid  
maintenance)



# Replicon: how the plasmid replicates

- A DNA sequence and associated factors
  - origin of replication, ori: a site on the DNA
  - protein and RNA factors that manage replication
- The replicon helps define the “ copy number” , the number of plasmid copies per cell

# A few common plasmid replicons

<u>PLASMID</u>	<u>REPLICON</u>	<u>COPY #</u>
pBR322	pMB1	15-20
pUC19	Modified form of pMB1 (RNA II mutation)	500-700
pACYC	p15A	18-22
pSC101	pSC101	~5

# Plasmid copy number is an important consideration


- High copy number plasmids
  - More copies/cell = easy to make a lot of DNA in a plasmid prep
  - Used for routine manipulation of small (<15 kb) recombinant DNAs
- Low copy number plasmids
  - For genes that are lethal or unstable in high copy number plasmids
  - Bacterial Artificial Chromosomes (BACs) that can propagate large (>100 kb) recombinant DNAs. Lots of DNA is a drag on cell physiology

# How to make sure a plasmid is maintained

- Plasmids contain **selectable markers**: genes carried by the plasmid that confer functions required for host survival
- Selection: only those cells with the plasmid will survive
  - Allows transformation (a rare event) to be easily detected
  - Cells cannot lose the plasmid, even if it causes a selective disadvantage (eg. slow growth or toxic gene product)

# Antibiotic resistance genes

- Antibiotic
  - **ampicillin** (related to penicillin, carbenicillin) inhibits cell wall synthesis
- Resistance gene
  - beta lactamase (bla) breaks down ampicillin, so cells carrying this gene are called **amp<sup>r</sup>**




## Keep in mind


- beta-lactamase is secreted into the medium where it breaks down the antibiotic and depletes it
- ampicillin sensitive cells can grow, defeating the selection

# Antibiotic resistance genes

- Antibiotic
  - chloramphenicol (cm) inhibits translation (peptidyl transferase activity of the ribosome)
- Resistance gene
  - chloramphenicol acetyl transferase (CAT) inactivates chloramphenicol




# Antibiotic resistance genes



- Antibiotic
  - tetracycline inhibits translation (interacts with ribosome 30S subunit)
- Resistance gene
  - Tet A (C ) protein confers resistance to by moving this antibiotic out of the cell


# Antibiotic resistance genes

- Antibiotic
  - **kanamycin**: aminoglycoside antibiotic, inhibits translation (interacts with ribosome 30S subunit)
  - others of this class include G418, neomycin
- Resistance gene
  - bacterial aminophosphotransferase inactivates kanamycin by phosphorylation





# pBR322: the first widely used plasmid



Transposition and restriction/ligation reactions were used to make it


## Using pBR322

Clone into sites in the Tc<sup>r</sup> gene,

Plate on ampicillin, then 'replica plate' on tetracycline.

Recombinants are amp resistant, tet sensitive


But: pBR322 has **low copy number, large size**



## pUC plasmids: improved cloning vectors

- Reduced size (about 2000 bp), easier for cells to handle
- Multiple cloning site (MCS, also called “ poly-linker” ): unique sites for lots of different restriction enzymes
- Very high copy number (mutation in RNA II)
- A “ blue-white” screening tool for recombinants (“ alpha complementation” is disrupted by foreign DNA in the MCS)


# Identifying recombinant plasmid: 'alpha complementation' helps



DNA cloned in the MCS interrupts the lacZ gene alpha fragment


# Alpha complementation

- **Plasmid** has the N-terminus of beta galactosidase (alpha fragment)
- Host *E. coli* has the C-terminus of beta galactosidase (omega fragment)
- Both the N- and C-terminal fragments are needed for enzyme activity
- if beta galactosidase is present, X-gal is cleaved, producing a bright blue product



# Alpha complementation on a petri plate

- Plasmid has alpha fragment of beta galactosidase with an MCS, bacterial host has omega fragment
- Clone DNA into the MCS, no alpha fragment
- No alpha fragment, no functional B-galactosidase
- No B-gal, no blue color
- white colonies



pUC19  
transformation  
plate

# Many specialized plasmids are available

- For expression of specific genes
  - Plasmids are engineered to include“ tags” on expressed protein, to assist in purification or tracking
- For cloning of unstable or toxic genes
  - Low copy number
- For cloning gigantic pieces of DNA
  - Low or single copy number

# Want a new plasmid? Addgene.org, a plasmid repository

- <https://www.addgene.org/>

# Keeping track of your plasmid sequence and structure


- A plasmid editor (ApE)
  - Save plasmid sequence
  - Include annotations to the plasmid
  - <http://biologylabs.utah.edu/jorgensen/wayned/ape/>



# How to create a recombinant plasmid (general protocol)

- 1) Prepare the plasmid vector
  - a) Make the plasmid linear (PCR, or restriction enzyme)
  - b) Purify, quantify the DNA
- 2) Prepare the DNA to be cloned
  - a) Obtain the DNA
    - PCR
    - Small pieces of genomic DNA
  - b) Treat the ends of the DNA to make them clonable
    - Primers with modified 5' ends
    - Add sequence to 3' ends with terminal transferase
    - Ligate adaptor sequences to ends
- 3) Combine the plasmid and clone DNA
  - a) Standard cloning: ligase
  - b) Gibson Assembly: Exonuclease, DNA pol, ligase
  - c) In vivo assembly (IVA): the cell assembles pieces with overlapping ends

**Gibson DNA assembly:** make synthetic genes, pathways, or *entire genomes*.




# Gibson assembly example

Cloning vector: plasmid pUC19

Cloned DNA: Pf1831 (*Pyrococcus furiosus* histone protein)

```
ATGGGAGAATTGCCAATTGCCCCAGTTGACAGGCTTATAAGAAAGGCAGG
TGCTGAAAGAGTTAGCGAGCAAGCAGCCAAAGTCCTAGCAGAGTACCTCG
AAGAGTACGCTATTGAGGTCGCAAAGAAGGCAGTAGAGTTTCGCAAGGCAC
GCAGGTAGAAAGACTGTTAAGGTTGAAGACATTAAGCTCGCAATTAAGAG
CTGA
```

pUC19 MCS



Cloning strategy: <http://nebuilder.neb.com/>


# Parameters for cloning by Gibson assembly:

The final product has to be circular

Need at least a 15 bp overlap at the cloning junctions (this is accomplished by extra sequences at the 5' end of the PCR primers)

Picking primers: can be done manually, or with an automated online tool

Overlaps	Oligo (Uppercase = gene-specific primer)	Anneals	F/R	3' Tm	3' Ta *
pUC19	tcgagctcggtacccGGCAATCTATTTGGAAATC	Pf1831	Fwd	56.0°C	56.0°C
pUC19	ttgcatgcctgcaggGATATACTTTAATTTCTGGGAGG	Pf1831	Rev	57.2°C	56.0°C




## "IVA" cloning: In Vivo Assembly

- *E. coli* has a recombination machinery that only requires short regions of homology (15-20 base pairs)
- Transform cells with DNA with overlaps, the recombination machinery will put them together automatically, and with reasonably high efficiency

**b**

**Sub-cloning**



You do

Cell does

# Moving DNA: transformation of *E.coli* with plasmid DNA

Cells acquire DNA by

- Chemical competence
- Electroporation
- Natural competence (in some lucky cases)

*E.coli* strain: must be antibiotic sensitive, and should not  
have restriction enzyme systems

Following transformation, the cells are cultured on a petri  
plate

## Making chemically competent *E. coli*

- Grow cells to “ mid-log” phase ( $A_{600}$  of 0.4) spin to get cell pellet
- Suspend cells in  $\text{CaCl}_2$  (100 mM), pellet again
- Suspend in small volume of  $\text{CaCl}_2$ /glycerol
- Freeze cells ( $-80^\circ\text{C}$ ) or go straight to transformation protocol

# Transformation of chemically competent cells

DNA binds to cells



- Mix DNA and competent cells, on ice for 30 min.

DNA uptake by cells



- Heat shock (42°C) for 1.5 minutes
- Add growth media, 37°C for 1 hour

Cells recover



- Plate on growth medium plus selection (antibiotic) for the plasmid

Selection occurs



Efficiency:  $\sim 10^5 - 10^8$  antibiotic resistant colonies/microgram plasmid DNA



# Transformation by electroporation

- Grow cells to  $A_{600}$  of 0.4
- Centrifuge and resuspend in water + 10% glycerol (do this 4 times to reduce conductivity)
- Place cells with DNA in electrode-containing cuvette, deliver electrical pulse
- If there is arcing (sparks) transformation efficiency will be poor (uneven transfer of charge). To avoid this make sure the ion concentration is very low (less than 10 mM salt)
- $> 10^9$  transformants/microgram DNA (ideally)

## Setting up a transformation – treatments and controls

1. *No plasmid* (negative control, nothing should grow on this plate)
2. *Supercoiled plasmid* of a known concentration (to determine efficiency of competent cells)
3. *Plasmid DNA* backbone: without clone DNA (background transformants)
4. *Plasmid DNA with clone DNA* (desired products)

# Ideal results from a transformation

- 1) No DNA--No colonies
- 2) 2 nanograms ( $2 \times 10^{-3}$  micrograms) supercoiled plasmid DNA--500 colonies (efficiency of cells:  $2.5 \times 10^5$  transformants per microgram DNA)
- 3) Plasmid alone--small number of colonies
- 4) Plasmid + insert—lots more colonies than #3

## Identify recombinant plasmids

- Screen colonies/plasmids for cloned DNA by PCR
- Alpha complementation: white colonies represent presence of insert DNA blocking functional beta galactosidase

Confirm clones by sequencing

# When cloning a piece of DNA consider:

1) Choice of plasmid: cloning sites? antibiotic? replicon?

2) Adding DNA to plasmid: how will the reaction be set up to give the desired product?

3) Moving DNA by transformation: what strain of *E. coli* will you transform into? Which method for transformation?

4) Screening for successful ligation products (recombinant plasmid DNA): how will the recombinant plasmids be identified?

# Basic gene cloning: plasmids and transformation, cutting and pasting

- Plasmid biology
- Adding DNA to a plasmid
- Transformation of *E. coli*

## Vectors for *E.coli* part II

- I. Bacteriophage (bacterial viruses):  
lambda and M13
  
- II. Moving and storing large DNA  
molecules: PACs, and BACs

## Readings:

- 1) 28 *MC4 Lambda and M13*. Short introduction to these historically important molecular biology tools, also two M13 protocols.
- 2) 33 *MC4 High capacity vectors*. Summary of vector types, with an emphasis on bacterial artificial chromosomes (BACs)
- 3) 29 *MC4 ~~Cre/Lox~~* Discussion of the cre-lox system for forcing site-specific recombination



# Bacteriophages: viral vectors for molecular cloning


- I. M13: a filamentous phage
  - Life cycle
  - genome structure
  
- II. Lambda: a “head and tail” phage
  - Life cycle
  - Basic cloning in lambda

# Bacteriophages


Viruses that infect bacteria

Morphologies:

- a) Head and tail
- b) Filamentous



(a) Head-and-tail



(b) Filamentous

- Nucleic acid molecule
  - Carrying genes for infection and replication
  - Surrounded by a protective protein coat (capsid)
- Infection (instead of transformation):
  - Phage attaches to outside of bacterium, injects DNA
  - Phage DNA is replicated, capsid proteins are made
  - Phage assembled and released

# Use of bacteriophages in cloning:

M13 -- single-stranded DNA genome

- Then:
  - Sequencing
  - Site-directed mutagenesis
- Now:
  - Phage display of foreign peptides


Lambda -- large-ish DNA fragments (25 kilobases)

- For gene cloning (large eukaryotic genes)
- Good selection schemes for recombinants
- Clone lots of precisely-sized DNA fragments for library construction

# M13: a filamentous bacteriophage

- Single-stranded, circular genome, 6.4 kb
- Can carry foreign DNA up to 6X the M13 genome size (36 kb) -- the larger the DNA, the less stable the clone
- Drawback: foreign DNA can be unstable (slow host cell growth – deletions confer a selective advantage)

M13 doesn't lyse cells, but it slows them down



"lawn"  
of *E. coli*

M13 infections form 'turbid' plaques

SS (+) strand

Infecting single-stranded circular viral DNA is converted to double-stranded replicating form (RF) DNA by host-encoded enzymes.

RF

ds

Isolate for cloning

Several rounds of replication occur through  $\theta$  structures.

The (-) strand of the RF DNA is transcribed into viral mRNAs.

The viral gene II product introduces a nick at a specific site (red arrow) in the (+) strand of the RF DNA.

Progeny (+) strand is synthesized continuously by movement of the replication machinery around the (-) strand template (rolling circle replication).

The completed progeny (+) strand is cleaved from the rolling circular structure by the viral gene II product (red arrow). The progeny strand then circularizes.


Synthesis of progeny (+) strand continues.

rolling circle structure


SS

progeny (+) strand

## M13 life cycle: an overview




# Engineered M13: a tool for cloning



alpha complementation (like pUC plasmids)



# M13 phage particle structure





# Bacteriophage lambda

- “ head and tail” phage, well-studied model
- Large, linear genome--48.5 kb
  - Central region of genome (“ stuffer” ) is dispensable for infectious growth--it can be engineered out
- Two infectious phases:
  - Lytic: replicating phase (killing/lysing cell)
  - Lysogenic: latent phase (integrated, waiting for opportunity)
- Can hold 5-25 kb DNA fragments

# Lambda as a cloning vector


## Phage capsid:

- Size of head is invariant
- Capsid must be filled
  - Too little DNA and it's not infectious
  - Too much DNA and the genome won't fit

a) Insertional vectors: can only increase genome size by 5% (about 3 to 5 kb)

b) Replacement vectors: remove “ stuffer” , can clone larger pieces of DNA, 8 to 24 kb (sufficient for many eukaryotic genes)

# Cloning in lambda phage – the genome view



- 1) Cut out stuffer, save right and left arms
- 2) Ligate with foreign DNA
- 3) “ Package” ligation mixture into phage heads
- 4) Plate mixture on *E. coli*, individual plaques represent recombinant clones

# Selection for *recombinant* lambda

There is a minimal size of DNA that can be packaged in lambda phage heads

- If the stuffer is removed the ligated “ arms” cannot be packaged (too small)
- Presence of recombinant DNA makes the phage large enough, so the recombinants are selected

# Cloning large DNA fragments

1) Make genomic *libraries*: the larger the DNA fragment, the fewer you need to make a complete *library* (representing the entire genome in fragments)

2) Clone DNA large enough to contain an entire eukaryotic gene

To get a functional gene, you need the whole gene and its regulatory regions

- Average human gene (exons and introns): 27 kilobases
- Regulatory regions may add another 10 to 100 kilobases upstream/downstream of gene

# Vectors for large recombinant DNA fragments

- Bacteriophage P1 plasmid: 70-100 kb
- YAC: 250-400 kb (or higher)
- PAC: 130-150 kb
- BAC: 120-300 kb (up to 700 kb is possible)


TABLE 1. High-capacity vectors for genomic cloning

Vector	Capacity (kb)	Replicon	Host	Copy number	Recovery of cloned DNA
P1	70–100	P1	<i>Escherichia coli</i>	1 (amplifiable)	Alkaline extraction
YAC	250–400	ARS	Yeast	1	Pulse-field gels
PAC	130–150	P1	<i>E. coli</i>	1	Alkaline extraction
BAC	120–300	F	<i>E. coli</i>	1	Alkaline extraction


BACs are most commonly used

## Phage P1 vectors:

clone large DNA fragments (85-100 kb)




# Phage P1 vectors



Efficiency of packaging is typically low: thus it is not good for making large genomic libraries



# Phage P1 vectors




*Cre* recombinase protein circularizes injected DNA at the *loxP* sites. DNA replicates using plasmid replicon. Plasmid copy number is increased by induction of P1 lytic replicon.


PACs: like P1 vectors but the DNA is not packaged (transfer by electroporation)

# Cre and Lox: a site-specific recombination system



The loxP site




Two 13 bp inverted repeats separated by an asymmetric linker



Binding of Cre protein to one or both *loxP* sites





Alignment of *loxP* sites



Asymmetric cleavage  
of *loxP* sites followed by  
strand invasion and synapse




kjihgfedcba



# BAC: Bacterial Artificial Chromosome

- F factor of E. coli:
  - 100 kb plasmid, propagates through conjugation
  - low copy number (1-2 copies per cell)
  - 2 genes (parA and parB): accurate partitioning during cell division
- BACs: engineered to have par genes, replication ori, cloning sites, selectable marker
- Holds very large pieces of DNA: up to 300 kb
- Fairly easy to manipulate: move into cells by transformation (electroporation)

# BAC vector map



## Vectors for *E.coli* part II

- I. Bacteriophage lambda and M13
- II. Moving and storing large DNA molecules: PACs and BACs